

# Veratox

## For DON HS (high sensitivity)




### Intended Use

Veratox for DON HS is a competitive direct ELISA that provides a quantitative analysis of lower levels of DON in cereal grains and grain products, such as wheat, wheat flour, wheat midds, wheat bran, corn, cornmeal, corn screenings, barley, malted barley, oats and processed cereal.

### Product Specifications

Lower limit of detection:	5 ppb
Range of quantitation:	25–250 ppb
Controls provided:	0, 25, 50, 100 and 250 ppb
Testing time:	20 minutes
Antibody cross-reactivity:	DON 100%. Also cross-reactive with 3-acetyl DON but does not bind to other trichothecenes.
Tests per kit:	Up to 38

### Results



### Ordering Information

Prod.#	Product description
8332	Veratox for DON HS
8052	Mycotoxin Extraction Kit

### Neogen Corporation

800/234-5333 • 517/372-9200 • E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com)

### Neogen Europe

+44 (0) 1292 525600 • E-mail: [info@neogeneurope.com](mailto:info@neogeneurope.com)